

aero

AERO offers a complete line of convenient wipe products that perform a wide variety of tasks.

CleanSwipe Wipes

DISINFECTING/DEODORIZING/CLEANING WIPES VIRUCIDAL*/BACTERICIDAL/TUBERCULOCIDAL/FUNGICIDAL/PSEUDOMONICIDAL For Use On Environmental Surfaces

Kills Pandemic 2009 H1N1 Influenza A virus (formerly called swine flu)

A non-woven disposable cloth containing a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces where control of cross-contamination is required.

This product is effective against the following pathogens in 2 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: Mycobacterium tuberculosis, *Human Hepatitis B virus, *Human Hepatitis C virus, *Herpes simplex virus 2 (genital Herpes virus), *HIV-1 (AIDS virus)**(At an exposure time of 1 min.), **See container label for complete HIV-1 (AIDS virus) information, *Influenza A (A2 Japan), Staphylococcus aureus, Pseudomonas aeruginosa, Salmonella enterica, Methicillin Resistant Staphylococcus aureus (MRSA), Escherichia coli O157:H7, Vancomycin Resistant Enterococcus faecalis (VRE).

This product is effective against the following pathogens in 5 minutes at 20°C/68°F, on hard, non-porous surfaces when used as directed: Trichophyton mentagrophytes (Athlete's foot fungus).

Areas of use: Hospitals, intensive care units, surgery, recovery, anesthesia, X-ray, Cat. lab., orthopedics, newborn nursery, respiratory therapy, emergency medical settings, day care, schools, school buses, laboratories, clinics, nursing facilities, medical offices, dental offices, veterinary facilities, janitorial, commercial, industrial facilities, spas and/or health clubs.

DIRECTIONS FOR USE: It is a violation of U.S. Federal Law to use this product in a manner inconsistent with its labeling.

This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes, but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high level disinfection. This product is not for use on the inhalation therapy and/or pulmonary diagnostic equipment.

DISINFECTION: Thoroughly wet pre-cleaned, hard, non-porous surface with a wipe, keep wet for 2 minutes (5 minutes if fungus is suspected), and allow to air dry. Use as many wipes as needed for the treated surface to remain wet for the entire contact time. For use as a general disinfectant wipe on hard, nonporous, inanimate surfaces such as stainless steel, Formica, glass tables, chairs, baskets, counter tops, exam tables, cages, toilets, basins, tubs, bed springs, frames, bedside, office furniture, chart racks, equipment, cabinets, door knobs, office machinery, sinks, mirrors, light fixtures, switches, carts, telephones, waste cans, laundry hampers, wheel chairs, walkers, urinals, plastic and other similar hard surfaces. You may use a wipe to pre-clean the surface to be disinfected, but thorough rewiping with a fresh wipe is required for disinfection. Use on food contact surfaces, must be followed by a potable water rinse.

DEODORIZING: Wipe hard, non-porous surface thoroughly to eliminate odors from bathrooms, fecal matter, urine, burns, smoke, mildew and mold, stale air, tobacco, sickrooms and musty closets.

CLEANING: To clean surfaces, wipe thoroughly with a fresh wipe, which may be used to deodorize at the same time as cleaning. Use in bathrooms, industrial plants, kitchens, sickrooms, musty closets and similar areas. When used to clean a hard non-porous surface prior to disinfection, a second wipe must be used for disinfection.

HIGHLIGHTS

- Convenient Pre-saturated Wipes
- Effective against- HIV-1, Hepatitis B and C, MRSA, VRE, Escherichia Coli, Pseudomonas aeruginosa, Salmonella enterica
- Kills Pandemic 2009 H1N1 influenza A virus (formerly called swine flu)

TYPICAL PROPERTIES

• ACTIVE INGREDIENTS:

N-Alkyl (60% C ₁₄ , 30% C ₁₆ , 5% C ₁₂ , 5% C ₁₈) Dimethyl Benzyl Ammonium Chloride.....	0.12%
N-Alkyl (68% C ₁₂ , 32% C ₁₄) Dimethyl Ethyl Benzyl Ammonium Chloride ..	0.12%
Isopropyl Alcohol.....	41.58%
OTHER INGREDIENTS†.....	58.18%

†Does not include weight of dry wipe.

- Appearance: 10 X 7 white cloth
- Fragrance: Alcohol
- Net Weight: 14 oz.

SAFETY CAUTIONS

Do not use as a baby wipe. Keep out of reach of children.

Read label and SDS before use.

Causes substantial, but temporary eye injury. Do not get in eyes or on clothing. Wear (specific appropriate protective eyewear such as goggles, face shield or safety glasses). Wash thoroughly with soap and water after handling and before eating, drinking or chewing gum, or using tobacco. Remove and wash contaminated clothing before reuse.